

High Speed Mini Centrifuge



Brand: HWLAB

Product Code: AC015HW

Availability: 2 - 3 Days

Call for Price: +2 01000 6546

Short Description

HWLAB® HW12 High speedmini Centrifuge is widely Used in physical and chemistry analysis...

Description

HWLAB® HW12 High speedmini Centrifuge is widely Used in physical and chemistry analysis, biochemistry,


Cellular and molecular biology, clinical Labs and blood donation centers.

features

Product features

- LCD display of time and speed
- 12x1.5ml/2.0ml capacity rotor
- Spin up to 15100xg (15000r/min)
- Metal rotor is autoclavable at 121°C in 20 minutes
- Adjustable speed from 800r/min to 15000r/min
- Double voltage, wide frequency, AC 110V/240V , 50-60Hz
- Separate Short Spin button for fast, convenient quick spins
- Excellent air flow minimizes heating-protects temperature-sensitive samples

- Lid automatically open at the end of the run to prevent sample warming and to reduce handing time
- Maintenance-free drive



Technical specifications

Technical specifications

Maximum rcf
Power supply
Rotor capacity
Speed
Timer

15100xg
 AC110V-240V/50-60Hz
 12x1.5ml/2.0ml, 2x8x0.2ml
 From 800r/min to 15000r/min
 15s to 99min

Dimension

23x22x13 cm

Dimension per carton

58x34x28 cm

Quantity per carton

2pcs

Weight

12KG

Product Gallery



R-12-16

