

ID-Card centrifuge



Brand: HWLAB

Product Code: AC014HW

Availability: 2 - 3 Days

Call for Price: +2 01000 6546

Short Description

HWLAB® ID-Card Centrifuge is widely Used in hospital,Clinical Labs...


Description

HWLAB® ID-Card Centrifuge is widely Used in hospital,Clinical Labs and blood donation centers. This is a specially Centrifugedeveloped for ID-System with a capacity of 24 ID-Cards.

features

Product features

- Maintenance-free brushless motor
 - LCD display of time and speed
 - Fixed speed 870r/min, rcf is 85xg
- Timer, from 15s to 99min, or continuous operation
 - Suitable for 24 ID-Cards.
- Humanized design. Three operating modes can choose ,including short spin/setting time/continuous operation
- Lid automatically open at the end of the run to prevent sample warming and to reduce handing time



Technical specifications

Technical specifications

Maximum rcf
Power supply
Rotor capacity

85xg
 AC100V-240V/50-60Hz
 24 ID-cards

Speed

870r/min

Timer

15s-99min / continuous operation

Motor type

Brushless motor

Dimension

370x330x175mm

Weight

11KG

Product Gallery

