

# Mini Clinical Centrifuge



**Brand:** HWLAB

**Product Code:** AC011HW

**Availability:** 2 - 3 Days

**Call for Price:** +2 01000 6546

## Short Description

HWLAB® mini Clinical Centrifuge is widely Used in physical and chemistry analysis


## Description

HWLAB® mini Clinical Centrifuge is widely Used in physical and chemistry analysis, biochemistry, Cellular and molecular biology, clinical Labs and blood donation centers.

## features

### Product features

- Maintenance-free brush less motor
- LCD display of time and speed
- Adjustable speed form 300r/min to 5000r/min
- Timer, form 15s to 99min, or continuous operation
- For 6x15ml/10ml/7ml/5ml tubes
- Lid automatically open at the end of the run to prevent sample warming and to reduce handing time



## Technical specifications

## Technical specifications

**Maximum rcf**  
**Power supply**  
**Rotor capacity**  
**Speed**  
**Timer**  
**Motor type**  
**Dimension**

2500xg  
 AC100V-240V/50-60Hz  
 6x15ml/10ml/7ml/5ml tubes  
 For 300r/min to 5000r/min  
 15s-99min / continuous operation  
 Brushless motor  
 212x270x158mm

**Weight**

4.5KG

**Product Gallery**

