

pipet4u® Revolution Electronic Pipette

Brand: AHN Biotechnologie GmbH

Product Code: AC037AH

Availability: 2 - 3 Days

Call for Price: +2 01000 6546



Short Description

Ultra bright OLED display for extra clear visibility...

Description

Comfortable and effortless pipetting

Ultra bright OLED display for extra clear visibility / Ergonomic grip shape / Ergonomic design for easy operation / Ergonomic pipetting button for easy operation setting and low strain pipetting (minimizes RSI)

Easy handling with one button - Control of all operating programs with the patented 1 button system

6 different operating programs:

Pipetting Mode

Dispensing Mode

Reverse Pipetting Mode

Mixing Mode

Auto Calibration Mode

Speed Control

Accurate and reliable pipetting - SCS-System controls piston movement for highest precision and accuracy

Easy self-calibration of pipette - ESM-System allows user calibration of pipette, e.g. when using different sample viscosities

Compatibility



[Compatibility with AHN myTip Standard and Macro Tips](#)



[Compatibility with AHN myTip FT Filter Tips](#)



[Compatibility with AHN myTip LT Tips](#)



[Compatibility with AHN myTip LFT Filter Tips](#)

Product Gallery

