

# Clinical Centrifuge



**Brand:** HWLAB  
**Product Code:** AC012HW  
**Availability:** 2 - 3 Days  
**Call for Price:** +2 01000 6546

## Short Description

HWLAB® Clinical Centrifuge is widely Used in physical and chemistry analysis, biochemistry...


## Description

HWLAB® Clinical Centrifuge is widely Used in physical and chemistry analysis, biochemistry, Cellular and molecular biology, clinical Labs and blood donation centers.

## features

### Product features

- Maintenance-free brushless motor
- LCD display of time and speed
- Adjustable speed form300r/min to 4000r/min
- Timer, form 15s to 99min, or continuous operation
- For 12x10ml/7ml/5ml blood collection tubes or12x15ml centrifuge tube or 8x50ml centrifuge tube
- Lid automatically open at the end of the run to prevent sample warming and to reduce handing time
- prevent sample warming and to reduce handing time



**Technical specifications**  
**Technical specification**

<b>Maximum rcf</b>	2460xg
<b>Power supply</b>	AC100V-240V/50-60Hz
<b>Rotor capacity</b>	12x15ml/10ml/7ml/5ml tubes or 8x50ml
<b>Speed</b>	Form 300r/min to 4000r/min
<b>Timer</b>	15s-99min / continuous operation
<b>Motor type</b>	Brushless motor
<b>Dimension</b>	370x330x175mm

**Weight**

**11KG**

**Product Gallery**

