

Mini Clinical Centrifuge



Brand: HWLAB

Product Code: AC011HW

Availability: 2 - 3 Days

Call for Price: +2 01000 6546

Short Description

HWLAB® mini Clinical Centrifuge is widely Used in physical and chemistry analysis


Description

HWLAB® mini Clinical Centrifuge is widely Used in physical and chemistry analysis, biochemistry, Cellular and molecular biology, clinical Labs and blood donation centers.

features

Product features

- Maintenance-free brush less motor
- LCD display of time and speed
- Adjustable speed form 300r/min to 5000r/min
- Timer, form 15s to 99min, or continuous operation
- For 6x15ml/10ml/7ml/5ml tubes
- Lid automatically open at the end of the run to prevent sample warming and to reduce handing time



Technical specifications

Technical specifications

Maximum rcf	2500xg
Power supply	AC100V-240V/50-60Hz
Rotor capacity	6x15ml/10ml/7ml/5ml tubes
Speed	For 300r/min to 5000r/min
Timer	15s-99min / continuous operation
Motor type	Brushless motor
Dimension	212x270x158mm

Weight

4.5KG

Product Gallery

