

Micropipette G model



Brand: HWLAB
Product Code: AC007HW
Availability: 2 - 3 Days
Call for Price: +2 01000 6546

Short Description

HWLAB® pipette G model is designed and manufactured to provide you with a range of robust

Description

HWLAB® pipette G model is designed and manufactured to provide you with a range of robust, accurate and precise pipettes. HWLAB® pipette G model is a fixed volume, air – displacement pipette.

features

Product features

- 121°C sterilization for low- half part.
- Very low operating force.
- Color indicates pipette volume.
- Positioned for perfect ergonomics.



Control button

- Very low operating force
- Colour indicates pipette volume
- Positioned for perfect ergonomics

Quick connection clip

- Remove lower part easily

Tip cone

- Tight fit to the tip

Ejector

- Very low operating force
- Positioned for perfect ergonomics

Product Gallery



Control button

- Very low operating force
- Colour indicates pipette volume
- Positioned for perfect ergonomics

Quick connection clip

- Remove lower part easily

Tip cone

- Tight fit to the tip

Ejector

- Very low operating force
- Positioned for perfect ergonomics